


LOQTEQ® Proximale Humerusplatte 3.5 • Proximal Humerus Plate 3.5


IC 6933-15 Instrumente • Instruments
IC 6933-35 Schraubenbank • Screw rack



ARTIKEL	ARTICLE	ANZAHL QUANTITY	ART.-NR. ART.-NO.
Leersieb Instrumente, LOQTEQ® PH 3.5	Tray for instruments, LOQTEQ® PH 3.5	1	IC 6933-16
Deckel für Leersiebe, groß	Lid for trays, large	1	IC 2008-00
1 Messinstrument für Schrauben ø3.5-4.0, bis L 90mm	Depth gauge for screws ø3.5-4.0, up to L 90mm	1	IS 7904-20
2 Bohrer Schnellkupplung ø2.3, L 180, Wendel 50	Twist drill ø2.3, L 180, coil 50, quick coupling	2	IU 7423-18
3 Bohrer Schnellkupplung ø2.7, L 150, Wendel 50	Twist drill ø2.7, L 150, coil 50, quick coupling	2	IU 7427-15
4 Bohrer Schnellkupplung ø3.5, L 110, Wendel 50	Twist drill ø3.5, L 110, coil 50, quick coupling	1	IU 7435-00
5 Schraubendrehereinsatz Duo, T15, Schnellkupplung	Screwdriver Duo, T15, quick coupling	2	IU 7825-56
6 Messstab LOQTEQ® PH, grün, für K-Draht L 150	Direct measuring device LOQTEQ®, green, for K-wire L 150	1	IU 7915-10
7 Spannbohrbuchse LOQTEQ® 3.5, Kompression 1mm	Load drill guide LOQTEQ® 3.5, compression 1mm	1	IU 8166-01
8 Spannbohrbuchse LOQTEQ® 3.5, Kompression 2mm	Load drill guide LOQTEQ® 3.5, compression 2mm	1	IU 8166-02
9 Grundhalterung für Spannbohrbuchse LOQTEQ® 3.5	Basic Insert for load drill guide LOQTEQ® 3.5	1	IU 8166-05
10 Bohrbuchse für Gleitloch LOQTEQ® 3.5, I-ø 2.8, rot	Drill guide for gliding hole LOQTEQ® 3.5, I-ø 2.8, red	2	IU 8166-10
11 Reduzierhülse für K-Draht ø1.6, grün	Reduction sleeve for K-wire ø1.6, green	2	IU 8166-15
12 Bohrbuchse für Rundloch LOQTEQ® 3.5, I-ø 2.4, grün	Drill guide for round hole LOQTEQ® 3.5, I-ø 2.4, green	2	IU 8166-30
13 K-Draht mit Trokarspitze, ø1.6, L 150	K-wire with trocar point, ø1.6, L 150	10	NK 0016-15
14 Behälter für K-Drähte L 150	Caddy for K-wire L 150	1	IC 0006-15
15 Kortikalisschraube 3.5, T15, selbstschneidend, L 12-38	Cortical screw 3.5, T15, self-tapping, L 12-38	4	SK 3514-xx-2
16 LOQTEQ® Kortikalisschraube 3.5, T15, selbstschneidend, L 18-44	LOQTEQ® Cortical screw 3.5, T15, self-tapping, L 18-44	4	SK 3525-xx-2
17 LOQTEQ® Spongiosaschraube 3.8, T15, L 28-60	LOQTEQ® Cancellous screw 3.8, T15, L 28-60	4	SP 3825-xx-2

LOQTEQ® Proximale Humerusplatte 3.5 • Proximal Humerus Plate 3.5

IC 6933-25 Plattenset • Set of Plates



ARTIKEL	ARTICLE	ANZAHL QUANTITY	ART.-NR. ART.-NO.
Leersieb Implantate, LOQTEQ® PH 3.5	Tray for implants LOQTEQ® PH 3.5, empty	1	IC 6933-26
Deckel für Leersiebe, groß	Lid for trays, large	1	IC 2008-00
1 Handgriff für Schnellkupplung groß, durchbohrt	Large handle, cannulated, quick coupling	1	IU 7706-00
2 Handgriff für Schnellkupplung m. Drehmomentbegrenzer 2.0 Nm	Handle with quick coupling, with torque limiter 2.0 Nm	1	IU 7707-20
3 Schraubenpinzette, selbsthaltend	Screw forceps, self-holding	1	IU 8004-00
4 Doppelbohrbuchse, Bohrer ø2.7/3.5 mit gefederter Zentrierung	Double drill guide ø2.7/3.5, with spring aided centering	1	IU 8116-60
5 Zielblock LOQTEQ® Proximale Humerusplatte 3.5	Aiming device LOQTEQ® Proximal Humerus Plate 3.5	1	IU 8176-01
6 Befestigungsschraube Zielblock LOQTEQ® KFI T15	Fixing screw aiming device LOQTEQ® SFI T15	1	IU 8176-03
7 LOQTEQ® Prox. Humerusplatte 3.5, 4-Loch, L 92, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 4 holes, L 92, Titanium	1	PH 3510-04-2
8 LOQTEQ® Prox. Humerusplatte 3.5, 5-Loch, L 105, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 5 holes, L 105, Titanium	1	PH 3510-05-2
9 LOQTEQ® Prox. Humerusplatte 3.5, 6-Loch, L 118, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 6 holes, L 118, Titanium	1	PH 3510-06-2
10 LOQTEQ® Prox. Humerusplatte 3.5, 8-Loch, L 143, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 8 holes, L 143, Titanium	1	PH 3510-08-2
OPTIONAL			
11 LOQTEQ® Prox. Humerusplatte 3.5, 4-Loch, L 92, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 4 holes, L 92, Titanium	1	PH 3510-04-2
12 LOQTEQ® Prox. Humerusplatte 3.5, 10-Loch, L 169, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 10 holes, L 169, Titanium	1	PH 3510-10-2
13 LOQTEQ® Prox. Humerusplatte 3.5, 12-Loch, L 195, Titan	LOQTEQ® Prox. Humerus Plate 3.5, 12 holes, L 195, Titanium	1	PH 3510-12-2
14 AcroPlate®, 3-Loch, R	AC-Plate, 3 holes, R	1	PX 3520-01-2
15 AcroPlate®, 3-Loch, L	AC-Plate, 3 holes, L	1	PX 3520-02-2
16 AcroPlate®, 4-Loch, R	AC-Plate, 4 holes, R	1	PX 3521-01-2
17 AcroPlate®, 4-Loch, L	AC-Plate, 4 holes, L	1	PX 3521-02-2